



ALPP Rabbit mAb

Catalog No	YP-rAb-17821
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB,IHC,IF,ELISA
Gene Name	ALPP;PLAP
Protein Name	Alkaline phosphatase, placental type
Purification Process	Protein A
Specificity	Endogenous
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Source	Monoclonal, Rabbit,IgG
Dilution	IHC 1:200-1:1000; WB 1:2000-1:10000; IF 1:200-1:1000; ELISA 1:5000-1:20000; Note: For IHC, we suggest antigen retrieval with TE buffer pH 9.0
Concentration	0.5 mg/ml
Purity	≥90%
Storage Stability	-15 ° C to -25 ° C/1 year(Do not lower than -25 ° C)
Synonyms	PLAP;Alkaline phosphatase, placental type;ALPP ; PLAP ; Alkaline phosphatase ; placental type ; Alkaline phosphatase Regan isozyme ; Placental alkaline phosphatase 1 ; PLAP-1;ALPP
Observed Band	70kD
Calculated Molecular Weight	58kD
Cell Pathway	Membrane
Tissue Specificity	Detected in placenta (at protein level).
Function	Catalytic activity:A phosphate monoester + H(2)O = an alcohol + phosphate.,cofactor: Binds 1 magnesium ion.,cofactor: Binds 2 zinc ions.,miscellaneous: In most mammals there are four different isozymes: placental, placental-like, intestinal and tissue non-specific (liver/bone/kidney).,online information:Alkaline phosphatase entry,polymorphism:Placental ALP is highly polymorphic, there are at least three common alleles.,similarity:Belongs to the alkaline phosphatase family.,subunit:Homodimer.,
Background	The protein encoded by this gene is an alkaline phosphatase, a metalloenzyme that catalyzes the hydrolysis of phosphoric acid monoesters. It belongs to a multigene family composed of four alkaline phosphatase isoenzymes. The enzyme functions as a homodimer and has a catalytic site containing one




magnesium and two zinc ions, which are required for its enzymatic function. The protein is primarily expressed in placental and endometrial tissue; however, strong ectopic expression has been detected in ovarian adenocarcinoma, serous cystadenocarcinoma, and other ovarian cancer cells. [provided by RefSeq, Jan 2015],

matters needing attention


Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Various whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-ALPP antibody. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody. Lane 1: HeLa Lane 2: JAR Lane 3: Rat small intestine Predicted band size: 58kDa Observed band size: 70kDa



Human placenta was stained with anti-ALPP rabbit antibody

