



ZNF160 Rabbit mAb

Catalog No	YP-rAb-18490
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB,FC
Gene Name	ZNF160,ZNF160;Zinc finger protein 160 (Zinc finger protein HZF5) (Zinc finger protein Kr18) (HKr18);KIAA1611
Protein Name	Zinc finger protein 160 (Zinc finger protein HZF5) (Zinc finger protein Kr18) (HKr18)
Purification Process	Protein A
Specificity	Endogenous
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Source	Monoclonal, Rabbit,IgG
Dilution	WB 1:1000-1:5000; FC 1:100-1:500;
Concentration	0.5 mg/ml
Purity	≥90%
Storage Stability	-15 ° C to -25 ° C/1 year(Do not lower than -25 ° C)
Synonyms	ZNF160;Zinc finger protein 160 (Zinc finger protein HZF5) (Zinc finger protein Kr18) (HKr18);KIAA1611
Observed Band	90kD
Calculated Molecular Weight	90kD
Cell Pathway	SUBCELLULAR LOCATION: Nucleus {ECO:0000305}.
Tissue Specificity	TISSUE SPECIFICITY: Ubiquitous. {ECO:0000269 PubMed:11410164}.
Function	May be involved in transcriptional regulation.
Background	The protein encoded by this gene is a Kruppel-related zinc finger protein which is characterized by the presence of an N-terminal repressor domain, the Kruppel-associated box (KRAB). The KRAB domain is a potent repressor of transcription; thus this protein may function in transcription regulation. Multiple transcript variants have been found for this gene. [provided by RefSeq, Apr 2016]
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Western blot analysis of lysates from SW480 cell, primary antibody was diluted at 1:1000, 4 ° over night, Dylight 800 secondary antibody was diluted at 1:10000, 37 ° 1 hour.

