




MS4A14 Rabbit mAb

Catalog No	YP-rAb-18492
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	MS4A14,MS4A13,MS4A16
Protein Name	Membrane-spanning 4-domains subfamily A member 14 (Testis development protein NYD-SP21)
Purification Process	Protein A
Specificity	Endogenous
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Source	Monoclonal, Rabbit,IgG
Dilution	WB 1:1000-1:5000;
Concentration	0.5 mg/ml
Purity	≥90%
Storage Stability	-15 ° C to -25 ° C/1 year(Do not lower than -25 ° C)
Synonyms	MS4A14;Membrane-spanning 4-domains subfamily A member 14 (Testis development protein NYD-SP21);MS4A13 MS4A16;MS4A13;MS4A16
Observed Band	75kD
Calculated Molecular Weight	75kD
Cell Pathway	SUBCELLULAR LOCATION: Membrane {ECO:0000305}; Multi-pass membrane protein {ECO:0000305}.
Tissue Specificity	
Function	May be involved in signal transduction as a component of a multimeric receptor complex. {ECO:0000250}.
Background	Predicted to be involved in cell surface receptor signaling pathway. Predicted to be located in membrane. Predicted to be active in plasma membrane. [provided by Alliance of Genome Resources, Apr 2025]
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Western blot analysis of lysates from MCF7 cell, primary antibody was diluted at 1:1000, 4 ° over night, Dylight 800 secondary antibody was diluted at 1:10000, 37 ° 1hour.

