



BAT5 Rabbit mAb

Catalog No	YP-rAb-18495
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	ABHD16A,BAT5,G5,NG26,PP199
Protein Name	Phosphatidylserine lipase ABHD16A (EC 3.1.-.-) (Alpha/beta hydrolase domain-containing protein 16A) (Abhydrolase domain-containing protein 16A) (HLA-B-associated transcript 5) (hBAT5) (Monoacylglycerol lipase ABHD16A) (EC 3.1.1.23) (Protein G5)
Purification Process	Protein A
Specificity	Endogenous
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Source	Monoclonal, Rabbit,IgG
Dilution	WB 1:1000-1:5000;
Concentration	0.5 mg/ml
Purity	≥90%
Storage Stability	-15 ° C to -25 ° C/1 year(Do not lower than -25 ° C)
Synonyms	BAT5;Phosphatidylserine lipase ABHD16A (EC 3.1.-.-) (Alpha/beta hydrolase domain-containing protein 16A) (Abhydrolase domain-containing protein 16A) (HLA-B-associated transcript 5) (hBAT5) (Monoacylglycerol lipase ABHD16A) (EC 3.1.1.23) (Protein G5);BAT5 G5 NG26;ABHD16A;G5;NG26;PP199
Observed Band	61kD
Calculated Molecular Weight	61kD
Cell Pathway	SUBCELLULAR LOCATION: Membrane {ECO:0000250 UniProtKB:Q9Z1Q2}; Multi-pass membrane protein {ECO:0000255}.
Tissue Specificity	
Function	Phosphatidylserine (PS) lipase that mediates the hydrolysis of phosphatidylserine to generate lysophosphatidylserine (LPS) (By similarity). LPS constitutes a class of signaling lipids that regulates immunological and neurological processes (By similarity). Has no activity towards diacylglycerol, triacylglycerol or lysophosphatidylserine lipase (PubMed:25290914). Also has monoacylglycerol lipase activity, with preference for 1-(9Z,12Z-octadecadienoyl)-glycerol (1-LG) and 2-glyceryl-15-deoxy-Delta(12,14)-prostaglandin J2 (15d-PGJ(2)-G) (PubMed:25290914). {ECO:0000250 UniProtKB:Q9Z1Q2, ECO:0000269 PubMed:25290914}.



Background


A cluster of genes, BAT1-BAT5, has been localized in the vicinity of the genes for tumor necrosis factor alpha and tumor necrosis factor beta. These genes are all within the human major histocompatibility complex class III region. The protein encoded by this gene is thought to be involved in some aspects of immunity. Alternatively spliced transcript variants have been described. [provided by RefSeq, Apr 2010]

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Western blot analysis of lysates from Caco-2 cell, primary antibody was diluted at 1:1000, 4 ° over night, Dylight 800 secondary antibody was diluted at 1:10000, 37 ° 1 hour.

