



PCSK9 Rabbit mAb

Catalog No	YP-rAb-18515
Isotype	IgG
Reactivity	Human,Rat
Applications	WB,IF,IP,ELISA
Gene Name	PCSK9 NARC1 PSEC0052
Protein Name	PCSK9
Purification Process	Protein A
Specificity	Endogenous
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Source	Monoclonal, Rabbit,IgG
Dilution	WB 1:2000-1:10000; IF 1:200-1:1000; ELISA 1:5000-1:20000; IP 1:50-1:200;
Concentration	0.5 mg/ml
Purity	≥90%
Storage Stability	-15 ° C to -25 ° C/1 year(Do not lower than -25 ° C)
Synonyms	Proprotein convertase subtilisin/kexin type 9 ; Neural apoptosis-regulated convertase 1 ; NARC-1 ; Proprotein convertase 9 ; PC9 ; Subtilisin/kexin-like protease PC9 ;
Observed Band	65kD(mature form),80kD(precursor)
Calculated Molecular Weight	74kD
Cell Pathway	Cytoplasm. Secreted. Endosome. Lysosome. Cell surface. Endoplasmic reticulum. Golgi apparatus. Autocatalytic cleavage is required to transport it from the endoplasmic reticulum to the Golgi apparatus and for the secretion of the mature protein. Localizes to the endoplasmic reticulum in the absence of LDLR and colocalizes to the cell surface and to the endosomes/lysosomes in the presence of LDLR. The sorting to the cell surface and endosomes is required in order to fully promote LDLR degradation.
Tissue Specificity	Expressed in neuro-epithelioma, colon carcinoma, hepatic and pancreatic cell lines, and in Schwann cells.
Function	urogenital system development, kidney development, liver development, regulation of receptor recycling, negative regulation of receptor recycling, regulation of receptor internalization, positive regulation of receptor internalization,proteolysis, neutral lipid metabolic process, acylglycerol metabolic process, triglyceride metabolic process, phospholipid metabolic process, glycerol ether metabolic process, induction of apoptosis, vacuolar transport, lysosomal

杭州臻优品生物科技有限公司

热销产品:

蛋白、一抗、抗体对、ELISA试剂盒、生化试剂盒
CCK8试剂盒、QPCR检测试剂盒

检测服务:

ELISA检测及定制服务 | 生化检测 | PCR、QPCR检测 | WB检测
ICO-IP检测 | 切片 | 染色 | 免疫组化 | 免疫荧光 | 透射电镜全套
| 宏基因组、转录组、基因组、蛋白组、代谢组测序



关注官网



关注客服



transport,steroid metabolic process, cholesterol metabolic process, macromolecule catabolic process, cellular response to starvation, response to endogenous stimulus, response to hormone stimulus, regulation of catabolic process,response to extracellular stimulus, response to organic substance, regulation of receptor activity, positive regulation of macromolecule metabolic process, negative regulation of macromolecule metabolic process, regulation of cell death,positive regulation of cell death, regulation of lipoprotein particle clearance, negative regulation of lipoprotein particle clearance, regulation of low-density lipoprotein particle clearance, negative regulation of low-density lipoprotein particle clearance, induction of programmed cell death, sterol metabolic process, protein processing, protein autoproccessing, organic ether metabolic process, organophosphate metabolic process, regulation of endocytosis,protein catabolic process, neuron differentiation, regulation of cellular catabolic process, response to nutrient levels,cellular response to extracellular stimulus, cellular response to nutrient levels, regulation of cellular protein metabolic process, low-density lipoprotein receptor metabolic process, receptor catabolic process, low-density lipoprotein receptor catabolic process, regulation of low-density lipoprotein receptor catabolic process, response to insulin stimulus, cellular response to insulin stimulus, cellular response to hormone stimulus, cellular response to stress,lipoprotein metabolic process, regulation of protein catabolic process, homeostatic process, response to starvation,cholesterol homeostasis, regulation of apoptosis, positive regulation of apoptosis, regulation of programmed cell death, positive regulation of programmed cell death, negative regulation of catalytic activity, receptor metabolic process, response to peptide hormone stimulus, regulation of neuron apoptosis, positive regulation of neuron apoptosis, negative regulation of molecular function, cellular protein catabolic process, cellular macromolecule catabolic process, positive regulation of endocytosis, glycerolipid metabolic process, intracellular transport, regulation of receptor-mediated endocytosis, positive regulation of receptor-mediated endocytosis, chemical homeostasis, positive regulation of transport, positive regulation of cellular component organization, negative regulation of multicellular organismal process, protein maturation, protein maturation by peptide bond cleavage, lipid homeostasis, sterol homeostasis, regulation of vesicle-mediated transport,

Background

cofactor:Calcium.,disease:Defects in PCSK9 are the cause of familial hypercholesterolemia 3 (FH3) [MIM:603776]. FH3 inheritance is autosomal dominant.,enzyme regulation:Inhibited by EGTA.,function:May be implicated in the differentiation of cortical neurons and may play a role in cholesterol homeostasis.,PTM:The soluble zymogen undergoes autocatalytic intramolecular processing in the endoplasmic reticulum, resulting in the cleavage of its propeptide that remains associated with the secreted enzyme.,similarity:Belongs to the peptidase S8 family.,similarity:Contains 1 peptidase S8 domain.,subunit:The precursor protein but not the mature protein may form multimers.,tissue specificity:Expressed in neuro-epithelioma, colon carcinoma, hepatic and pancreatic cell lines, and in Schwann cells.,


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.





Various whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-PCSK9 antibody. The HRP-conjugated Goat anti-Rabbit IgG (H + L) antibody was used to detect the antibody. Lane 1: HeLa Lane 2: A431 Lane 3: Rat liver
Predicted band size: 74kDa Observed band size: 65(mature form),80(precu

