



NCKPL rabbit pAb

Catalog No	YP-Ab-09762
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	NCKAP1L HEM1
Protein Name	Nck-associated protein 1-like (Hematopoietic protein 1) (Membrane-associated protein HEM-1)
Immunogen	Synthesized peptide derived from human NCKPL AA range: 582-632
Specificity	This antibody detects endogenous levels of human NCKPL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	124kD
Cell Pathway	Cell membrane ; Single-pass membrane protein ; Cytoplasmic side . Cytoplasm . Localizes to the leading edge of polarized neutrophils. .
Tissue Specificity	Expressed only in cells of hematopoietic origin (PubMed:7643388, PubMed:1932118). Expressed in neutrophils (at protein level) (PubMed:16417406). Expressed in T-cells (at protein level) (PubMed:32647003).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images