



CFDP1 Rabbit pAb

Catalog No	YP-Ab-19438
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	CFDP1 BCNT
Immunogen	Synthesized peptide derived from human CFDP1. AA range:142-242
Specificity	This antibody detects endogenous levels of CFDP1 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Synonyms	CFDP1 BCNT;Craniofacial development protein 1 (Bucentaur)
Observed Band	
Calculated Molecular Weight	33kD
Cell Pathway	Chromosome , centromere , kinetochore .
Tissue Specificity	Ubiquitous.
Function	May play a role during embryogenesis.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

