



IGF-1 Rabbit pAb

Catalog No	YP-Ab-19441
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB、ELISA
Gene Name	IGF1 IBP1
Immunogen	Synthesized peptide derived from human IGF1. AA range:1-130
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat IGF1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:1000-2000; ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Synonyms	IGF1 IBP1;Insulin-like growth factor I;IBP1 ; IGF-IA ; IGF1A ; IGFI ; insulin-like growth factor 1 ; insulin-like growth factor IA ; mechano growth factor ; MGF ; somatomedin C
Observed Band	
Calculated Molecular Weight	21kD
Cell Pathway	Secreted .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

