



ARMC5 mouse mAb

Catalog No	YP-mAb-19442
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	ARMC5
Immunogen	Synthesized peptide derived from human ARMC5. AA range:533-603
Specificity	This antibody detects endogenous levels of ARMC5 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Synonyms	ARMC5;Armadillo repeat-containing protein 5
Observed Band	
Calculated Molecular Weight	103kD
Cell Pathway	Cytoplasm .
Tissue Specificity	
Function	Involved in fetal development , T-cell function and adrenal gland growth homeostasis (By similarity) . Negatively regulates adrenal cells survival. Plays a role in steroidogenesis , modulates steroidogenic enzymes expression and cortisol production .
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images