




PLEK mouse mAb

Catalog No	YP-mAb-11483
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	PLEK P47
Protein Name	PLEK
Immunogen	Synthesized peptide derived from human PLEK AA range: 99-149
Specificity	This antibody detects endogenous levels of PLEK at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	39kD
Cell Pathway	extracellular region,cytoplasm,cytosol,membrane,ruffle membrane,
Tissue Specificity	
Function	function:Major protein kinase C substrate of platelets, its exact function is not known.,similarity:Contains 1 DEP domain.,similarity:Contains 2 PH domains.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using PLEK mouse mAb