



MFRN1 rabbit pAb

Catalog No	YP-Ab-10173
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SLC25A37 MFRN MSCP HT015
Protein Name	Mitoferrin-1 (Mitochondrial iron transporter 1) (Mitochondrial solute carrier protein) (Solute carrier family 25 member 37)
Immunogen	Synthesized peptide derived from human MFRN1 AA range: 187-237
Specificity	This antibody detects endogenous levels of human MFRN1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	37kD
Cell Pathway	Mitochondrion inner membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

