




beta Tubulin Rabbit mAb

Catalog No	YP-nAb-1176
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB,IHC-P,ICC/IF,FC
Alternative Name	Beta 4 tubulin; TBB5; TUBB2; TUBB2A; tubulin beta 2A; beta tubulin; tubulin beta chain; Tubulin beta 2b; Alpha tubulin; Alpha-tubulin ubiquitous;
Antibody dilution ratio	WB-1:1000-1:5000; IHC-P-1:20-1:50; ICC/IF-1:200-1:500; FC-1:20
SwissProt ID	P07437
Host Species	Rabbit
Conjugation	Unconjugated
Modification	Unmodified
Gene ID	203068
Gene Name	TUBB
Clonality	Monoclonal Antibody
Molecular Weight	Calculated MW:50 kDa ;Observed MW: 55 kDa
Immunogen	A synthesized peptide derived from human beta Tubulin
Research Field	Isotype>Loading Controls
Background	Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable site on the beta chain and one at a non-exchangeable site on the alpha chain.
Purification Method	Affinity Purified
Form	Liquid
Storage Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.
Storage Conditions	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.







Products Images



Western blot analysis of beta Tubulin in K562, rat Brain, C6, 3T3 lysates using beta Tubulin antibody.



Immunocytochemistry analysis of beta Tubulin (green) in HeLa using beta Tubulin antibody, and DAPI (blue)



Immunohistochemistry analysis of paraffin-embedded Human brain using beta Tubulin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

杭州臻优品生物科技有限公司

热销产品:

蛋白、一抗、抗体对、ELISA试剂盒、生化试剂盒
CCK8试剂盒、QPCR检测试剂盒

检测服务:

ELISA检测及定制服务 | 生化检测 | PCR、QPCR检测 | WB检测
ICO-IP检测 | 切片 | 染色 | 免疫组化 | 免疫荧光 | 透射电镜全套
| 宏基因组、转录组、基因组、蛋白组、代谢组测序



关注官网



关注客服