




beta Tubulin (6D8) Mouse mAb

Catalog No	YP-nAb-1190
Isotype	IgG1
Reactivity	Human,Mouse,Rat,Hamster,Chicken,Dog,Monkey,Rabbit,Sheep,Insect,Yeast
Applications	WB,IHC-F,IHC-P,ICC/IF
Alternative Name	Beta 4 tubulin; TBB5; TUBB2; TUBB2A; tubulin beta 2A; beta tubulin; tubulin beta chain; Tubulin beta 2b; Alpha tubulin; Alpha-tubulin ubiquitous;
Antibody dilution ratio	WB-1:5000-1:10000; IHC-1:50-1:200; IF-1:100-1:200
SwissProt ID	P07437
Host Species	Mouse
Conjugation	Unconjugated
Modification	Unmodified
Gene ID	203068
Gene Name	TUBB
Clonality	Monoclonal Antibody
Molecular Weight	Calculated MW: 50 kDa; Observed MW: 55 kDa
Immunogen	Synthetic Peptide of β -Tubulin
Research Field	Isotype/Loading Controls
Background	Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable site on the beta chain and one at a non-exchangeable site on the alpha chain.
Purification Method	Affinity Purified
Form	Liquid
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Storage Conditions	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.






Products Images



Western blot analysis of betaTubulin in CHO-K1, Jurkat, 3T3 and C6 lysates using betaTubulin antibody.



Immunofluorescence analysis of beta Tubulin in Hela using beta Tubulin antibody(Left) and DAPI (Right).




Immunohistochemistry analysis of paraffin-embedded Human colon tissue using u3b2tubulin antibody(5G3).High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.





Western blot analysis of β Tubulin (6D8) in HeLa, A549, K562, MCF-7 lysates using β Tubulin (6D8) antibody.



杭州臻优品生物科技有限公司

热销产品:

蛋白、一抗、抗体对、ELISA试剂盒、生化试剂盒
CCK8试剂盒、QPCR检测试剂盒

检测服务:

ELISA检测及定制服务 | 生化检测 | PCR、QPCR检测 | WB检测
ICO-IP检测 | 切片 | 染色 | 免疫组化 | 免疫荧光 | 透射电镜全套
| 宏基因组、转录组、基因组、蛋白组、代谢组测序



关注官网



关注客服