




Annexin X Monoclonal Antibody

Catalog No	YP-mAb-03705
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ANXA10
Protein Name	Annexin A10
Immunogen	The antiserum was produced against synthesized peptide derived from human ANXA10. AA range:261-310
Specificity	Annexin X Monoclonal Antibody detects endogenous levels of Annexin X protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ANXA10; ANX14; Annexin A10; Annexin-10; Annexin-14
Observed Band	38kD
Cell Pathway	mitochondrion,
Tissue Specificity	Colon,Lung,
Function	similarity:Belongs to the annexin family.,similarity:Contains 4 annexin repeats.,
Background	This gene encodes a member of the annexin family. Members of this calcium-dependent phospholipid-binding protein family play a role in the regulation of cellular growth and in signal transduction pathways. The function of this gene has not yet been determined. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using Annexin X Monoclonal Antibody