



ABI3BP Mouse mAb

Catalog No	YP-mAb-19453
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	ABI3BP NESHBP TARSH
Immunogen	Synthesized peptide derived from human ABI3BP. AA range:146-216
Specificity	This antibody detects endogenous levels of ABI3BP at Human,Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Synonyms	ABI3BP NESHBP TARSH;Target of Nesh-SH3 (Tarsh) (ABI gene family member 3-binding protein) (Nesh-binding protein) (NeshBP)
Observed Band	
Calculated Molecular Weight	117kD
Cell Pathway	Secreted .
Tissue Specificity	Expressed in brain , heart , lung , liver , pancreas kidney and placenta.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

