



UBP2 rabbit pAb

Catalog No	YP-Ab-09622
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	USP2 UBPA1
Protein Name	Ubiquitin carboxyl-terminal hydrolase 2 (EC 3.4.19.12) (41 kDa ubiquitin-specific protease) (Deubiquitinating enzyme 2) (Ubiquitin thioesterase 2) (Ubiquitin-specific-processing protease 2)
Immunogen	Synthesized peptide derived from human UBPA AA range: 452-502
Specificity	This antibody detects endogenous levels of human UBPA
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	67kD
Cell Pathway	Cytoplasm . Cytoplasm, perinuclear region . Localizes in the spermatid head in late-elongating spermatids in the thin area between the outer acrosomal membrane and the plasma membrane. .; [Isoform 4]: Nucleus . Membrane ; Peripheral membrane protein . Cytoplasm . Predominantly expressed at membranes. .
Tissue Specificity	Expressed in mesangial cells of the kidney and in different types of glomerulonephritides (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images