



SULF2 rabbit pAb

Catalog No	YP-Ab-09758
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	SULF2 KIAA1247 UNQ559/PRO1120
Protein Name	Extracellular sulfatase Sulf-2 (hSulf-2) (EC 3.1.6.-)
Immunogen	Synthesized peptide derived from human SULF2 AA range: 690-740
Specificity	This antibody detects endogenous levels of human SULF2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	96kD
Cell Pathway	Endoplasmic reticulum . Golgi apparatus, Golgi stack . Cell surface . Also localized on the cell surface. .
Tissue Specificity	Expressed at highest levels in the ovary, skeletal muscle, stomach, brain, uterus, heart, kidney and placenta.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

