



EPHX2 mouse mAb

Catalog No	YP-mAb-17670
Isotype	IgG
Reactivity	Human, Mouse,Rat
Applications	WB
Gene Name	EPHX2
Protein Name	
Immunogen	
Specificity	EPHX2
Formulation	Bifunctional epoxide hydrolase 2 [Includes: Cytosolic epoxide hydrolase 2 (CEH) (EC 3.3.2.10) (Epoxide hydratase) (Soluble epoxide hydrolase) (SEH); Lipid-phosphate phosphatase (EC 3.1.3.76)] Human Gene Id: 2053 Human Swiss Prot No: P34913 Mouse Gene Id: 13850 Mouse Swiss Prot No: P34914 Rat Gene Id: 65030 Rat Swiss Prot No: P80299
Source	Mouse, Monoclonal
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	The antibody was affinity-purified from Mouse serum by affinity-chromatography using specific immunogen.



Observed Band 61kD

Cell Pathway

Tissue Specificity

Function

Background

matters needing attention Avoid repeated freezing and thawing!

Usage suggestions This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

