



WDR6 Rabbit pAb

Catalog No	YP-Ab-19460
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	WDR6
Immunogen	Synthesized peptide derived from human WDR6. AA range:889-959
Specificity	This antibody detects endogenous levels of WDR6 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Synonyms	WDR6;WD repeat-containing protein 6
Observed Band	
Calculated Molecular Weight	123kD
Cell Pathway	Cytoplasm . Colocalizes in the cytoplasm with STK11/LKB1.
Tissue Specificity	Ubiquitous.
Function	Enhances the STK11/LKB1-induced cell growth suppression activity. Negative regulator of amino acid starvation-induced autophagy.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images