




CD66b (PTR1477) Mouse mAb

Catalog No	YP-Ab-19463
Isotype	IgG
Reactivity	Human
Applications	IF、IHC
Gene Name	CEACAM8 CGM6
Immunogen	The specific immunogen used to produce this antibody is proprietary information.
Specificity	Endogenous
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Source	Monoclonal,Mouse,IgG
Purification	Recombinant Antibody expressed in animal component-free (ACF) media, purified via Protein A affinity chromatography.
Dilution	IHC 1:200-1:1000; IF 1:200-1:1000 Note: For IHC, we suggest antigen retrieval with TE buffer pH 9.0 (Cat#YS0004)
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Synonyms	CEACAM8 CGM6;Carcinoembryonic antigen-related cell adhesion molecule 8 (CD67 antigen) (Carcinoembryonic antigen CGM6) (Non-specific cross-reacting antigen NCA-95) (CD antigen CD66b);Carcinoembryonic antigen CGM6 antibody ; Carcinoembryonic antigen gene family member 6 antibody ; Carcinoembryonic antigen related cell adhesion molecule 8 antibody ; Carcinoembryonic antigen-related cell adhesion molecule 8 antibody ; CD 66b antibody ; CD 67 antibody ; CD66b antibody ; CD66b antigen antibody ; CD67 antibody ; CD67 antigen antibody ; CEACAM 8 antibody ; CEACAM8 antibody ; CEAM8_HUMAN antibody ; CGM 6 antibody ; CGM6 antibody ; NCA 95 antibody ; NCA95 antibody ; Non-specific cross-reacting antigen NCA-95 antibody ; Nonspecific cross reacting antigen NCA 95 antibody ; Nonspecific cross reacting antigen NCA95 antibody ; CD66b nanobody ;
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Lipid-anchor , GPI-anchor . Cell surface .
Tissue Specificity	Expressed in leukocytes of chronic myeloid Leukemia patients and bone marrow.



Function	similarity:Belongs to the immunoglobulin superfamily. CEA family. ,similarity:Contains 1 Ig-like V-type (immunoglobulin-like) domain. ,similarity:Contains 2 Ig-like C2-type (immunoglobulin-like) domains. ,tissue specificity:Expressed in leukocytes of chronic myeloid Leukemia patients and bone marrow. ,
Background	similarity:Belongs to the immunoglobulin superfamily. CEA family. ,similarity:Contains 1 Ig-like V-type (immunoglobulin-like) domain. ,similarity:Contains 2 Ig-like C2-type (immunoglobulin-like) domains. ,tissue specificity:Expressed in leukocytes of chronic myeloid Leukemia patients and bone marrow. ,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Human spleen was stained with anti-CD66b (PTR1477) Mouse antibody