




RNF180 Mouse mAb

Catalog No	YP-mAb-19462
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	
Immunogen	Recombinant protein (or fragment).This information is considered to be commercially sensitive.
Specificity	
Formulation	
Source	Monoclonal,Mouse,IgG
Purification	Affinity purification
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Synonyms	RINES; RNF180
Observed Band	65kDa
Calculated Molecular Weight	68kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	Predicted to enable ubiquitin conjugating enzyme binding activity and ubiquitin protein ligase activity. Predicted to be involved in norepinephrine metabolic process; positive regulation of proteasomal ubiquitin-dependent protein catabolic process; and serotonin metabolic process. Predicted to act upstream of or within several processes, including adult behavior; positive regulation of protein ubiquitination; and protein polyubiquitination. Predicted to be located in nuclear envelope. Predicted to be integral component of membrane. Predicted to be intrinsic component of endoplasmic reticulum membrane.
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of various lysates using RNF180 Mouse mAb (A15241) at 1:1000 dilution.

Secondary antibody: HRP-conjugated Goat anti-Mouse IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 3s.