



UN13A mouse mAb

Catalog No	YP-mAb-09279
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	UNC13A KIAA1032
Protein Name	Protein unc-13 homolog A (Munc13-1)
Immunogen	Synthesized peptide derived from human UN13A AA range: 35-85
Specificity	This antibody detects endogenous levels of human UN13A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	187kD
Cell Pathway	Cytoplasm . Cell membrane ; Peripheral membrane protein . Cell junction, synapse, presynaptic cell membrane ; Peripheral membrane protein . Cell junction, synapse, presynaptic active zone . Translocated to the plasma membrane in response to phorbol ester binding. .
Tissue Specificity	Expressed in pancreatic islet cells (PubMed:12871971). Expressed in melanocytes (PubMed:23999003).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images