




LIX1 Monoclonal Antibody

Catalog No	YP-mAb-07709
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LIX1 C5orf11
Protein Name	Protein limb expression 1 homolog
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	LIX1 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	31kD
Cell Pathway	cytoplasm,
Tissue Specificity	Brain,
Function	similarity:Belongs to the LIX1 family.,
Background	similarity:Belongs to the LIX1 family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using LIX1 Monoclonal Antibody