



GALT Rabbit pAb

Catalog No	YP-Ab-19472
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	GALT
Immunogen	Synthesized peptide derived from human GALT. AA range:103-203
Specificity	This antibody detects endogenous levels of GALT at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Synonyms	GALT;Galactose-1-phosphate uridylyltransferase (Gal-1-P uridylyltransferase) (UDP-glucose--hexose-1-phosphate uridylyltransferase)
Observed Band	
Calculated Molecular Weight	42kD
Cell Pathway	
Tissue Specificity	
Function	Plays an important role in galactose metabolism.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

