




OGFOD1 Mouse mAb

Catalog No	YP-mAb-19468
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	
Immunogen	Recombinant protein (or fragment).This information is considered to be commercially sensitive.
Specificity	
Formulation	
Source	Monoclonal,Mouse,IgG
Purification	Affinity purification
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Synonyms	TPA1; OGFOD1
Observed Band	63kDa
Calculated Molecular Weight	63kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	Enables peptidyl-proline 3-dioxygenase activity. Involved in several processes, including peptidyl-proline hydroxylation; regulation of translational termination; and stress granule assembly. Located in cytoplasmic stress granule; cytosol; and nucleoplasm.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HeLa cells, using OGFOD1 Mouse mAb (A16543) at 1:1000 dilution.
Secondary antibody: HRP-conjugated Goat anti-Mouse IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25μg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 5s.