



Willin Rabbit pAb

Catalog No	YP-Ab-19488
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	FRMD6 C14orf31
Immunogen	Synthesized peptide derived from human Willin. AA range:1-130
Specificity	This antibody detects endogenous levels of Willin at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Synonyms	FRMD6 C14orf31;FERM domain-containing protein 6 (Willin)
Observed Band	
Calculated Molecular Weight	68kD
Cell Pathway	Cytoplasm . Cell membrane ; Peripheral membrane protein ; Cytoplasmic side . Can colocalize with actin.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images