



DMKN Mouse mAb

Catalog No	YP-mAb-19483
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	DMKN UNQ729/PRO1411
Immunogen	Synthesized peptide derived from human DMKN. AA range:298-398
Specificity	This antibody detects endogenous levels of DMKN at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Synonyms	DMKN UNQ729/PRO1411;Dermokine (Epidermis-specific secreted protein SK30/SK89)
Observed Band	
Calculated Molecular Weight	52kD
Cell Pathway	Secreted .
Tissue Specificity	Expressed in epidermis; in the spinous and granular layers and in placenta. Also found in the epithelia of the small intestine , macrophages of the lung and endothelial cells of the lung. Isoform 15 is expressed in epidermis and placenta. Isoform 1 is expressed in epidermis.
Function	May act as a soluble regulator of keratinocyte differentiation.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images