



TPPP3 Rabbit pAb

Catalog No	YP-Ab-19491
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	TPPP3 CGI-38
Immunogen	Synthesized peptide derived from human TPPP3. AA range:17-117
Specificity	This antibody detects endogenous levels of TPPP3 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Synonyms	TPPP3 CGI-38;Tubulin polymerization-promoting protein family member 3 (TPPP/p20);P20
Observed Band	
Calculated Molecular Weight	19kD
Cell Pathway	Cytoplasm . Cytoplasm , cytoskeleton .
Tissue Specificity	Expressed in endometrium during the mid-secretory phase (LH + 7) (at protein level) .
Function	Regulator of microtubule dynamic that has microtubule bundling activity . Required for embryo implantation; possibly by regulating beta-catenin (By similarity) . Also required for decidualization via regulation of beta-catenin .
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images