



Myosin Heavy Chain, Smooth Muscle (SMMHC) (PT1919R) Rabbit mAb

Catalog No	YP-rAb-18546
Isotype	IgG
Reactivity	Human, Mouse, Rat
Applications	WB、IHC、IF、ELISA
Gene Name	MYH11 KIAA0866
Immunogen	The specific immunogen used to produce this antibody is proprietary information.
Specificity	Endogenous
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% BSA
Source	Monoclonal, Rabbit, IgG
Purification	Recombinant Antibody expressed in animal component-free (ACF) media, purified via Protein A affinity chromatography.
Dilution	IHC 1:200-1:1000; WB 1:2000-1:10000; IF 1:200-1:1000; ELISA 1:5000-1:20000 Note: For IHC, we suggest antigen retrieval with TE buffer pH 9.0 (Cat#YS0004)
Concentration	0.5 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year (Do not lower than -25°C)
Synonyms	MYH11 KIAA0866; Myosin Heavy Chain, Smooth Muscle; Myosin-11 ; Myosin heavy chain 11 ; Myosin heavy chain, smooth muscle isoform ; SMMHC ;
Observed Band	255kD
Calculated Molecular Weight	227kD
Cell Pathway	Cytoplasmic
Tissue Specificity	Smooth muscle; expressed in the umbilical artery , bladder , esophagus and trachea. Isoform 1 is mostly found in slowly contracting tonic muscles.
Function	Disease: A chromosomal aberration involving MYH11 is found in acute myeloid leukemia of M4EO subtype. Pericentric inversion inv (16) (p13;q22) . The inversion produces a fusion protein consisting of the 165 N-terminal residues of CBF-beta (PEPB2) and the tail region of MYH11. , Disease: Defects in MYH11 are the cause of aortic aneurysm familial thoracic type 4 (AAT4) [MIM:132900]; also known as familial thoracic aortic aneurysm and dissection (TAAD) . Aneurysms and dissections of the aorta usually result from degenerative changes in the aortic wall. Thoracic aortic aneurysms and dissections are primarily associated with a characteristic histologic appearance known as 'medial necrosis' or 'Erdheim cystic medial necrosis' in which there is degeneration and fragmentation of elastic fibers , loss of smooth muscle cells , and an accumulation of basophilic ground substance. Patients with AAT4 show marked aortic stiffness. Pathological aortas



show large areas of medial degeneration with very low smooth muscle cells content. ,Domain:The rodlike tail sequence is highly repetitive , showing cycles of a 28-residue repeat pattern composed of 4 heptapeptides , characteristic for alpha-helical coiled coils. ,Function:Muscle contraction. ,miscellaneous:Each myosin heavy chain can be split into 1 light meromyosin (LMM) and 1 heavy meromyosin (HMM) . It can later be split further into 2 globular subfragments (S1) and 1 rod-shaped subfragment (S2) . ,similarity:Contains 1 IQ domain. ,similarity:Contains 1 myosin head-like domain. ,subcellular location:Identified by mass spectrometry in melanosome fractions from stage I to stage IV. Thick filaments of the myofibrils. ,subunit:Muscle myosin is a hexameric protein that consists of 2 heavy chain subunits (MHC) , 2 alkali light chain subunits (MLC) and 2 regulatory light chain subunits (MLC-2) . ,tissue specificity:Smooth muscle; expressed in the umbilical artery , bladder , esophagus and trachea. ,

Background

The protein encoded by this gene is a smooth muscle myosin belonging to the myosin heavy chain family. The gene product is a subunit of a hexameric protein that consists of two heavy chain subunits and two pairs of non-identical light chain subunits. It functions as a major contractile protein , converting chemical energy into mechanical energy through the hydrolysis of ATP. The gene encoding a human ortholog of rat NUDE1 is transcribed from the reverse strand of this gene , and its 3' end overlaps with that of the latter. The pericentric inversion of chromosome 16 [inv (16) (p13q22)] produces a chimeric transcript that encodes a protein consisting of the first 165 residues from the N terminus of core-binding factor beta in a fusion with the C-terminal portion of the smooth muscle myosin heavy chain. This chromosomal rearrangement is associated with acute myeloid leukemia of the M4Eo subtype. Alter


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images




Various whole cell lysates were separated by 4-8% SDS-PAGE, and the membrane was blotted with anti-Myosin Heavy Chain, Smooth Muscle (SMMHC) (PT1919R) antibody. The HRP-conjugated Goat anti-Rabbit IgG (H + L) antibody was used to detect the antibody.


Lane 1: Rat spleen

Predicted band size: 227kDa


Observed band size: 255kDa




Mouse colon was stained with anti-Myosin Heavy Chain, Smooth Muscle (SMMHC) (PT1919R) Rabbit antibody



Rat colon was stained with anti-Myosin Heavy Chain, Smooth Muscle (SMMHC) (PT1919R) Rabbit antibody



Human colon was stained with anti-Myosin Heavy Chain, Smooth Muscle (SMMHC) (PT1919R) Rabbit antibody



Human appendix was stained with anti-Myosin Heavy Chain, Smooth Muscle (SMMHC) (PT1919R) Rabbit antibody