



NBPF5 Monoclonal Antibody

Catalog No	YP-mAb-02692
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NBPF5
Immunogen	The antiserum was produced against synthesized peptide derived from human NBPF5. AA range:302-351
Specificity	NBPF5 Monoclonal Antibody detects endogenous levels of NBPF5 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NBPF5;Neuroblastoma breakpoint family member 5
Observed Band	41kD
Calculated Molecular Weight	41kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Expressed in brain and medulla.
Function	
Background	NBPF5 (neuroblastoma breakpoint family member 5) is a 351 amino acid cytoplasmic protein that is expressed in medulla and brain and belongs to the NBPF family. NBPF5 contains one NBPF domain and is encoded by a gene that maps to human chromosome 1p13. Chromosome 1 spans 260 million base pairs, contains over 3000 genes and comprises nearly 8% of the human genome. Chromosome 1 houses a large number of disease-associated genes, including those that are involved in familial adenomatous polyposis, Stickler syndrome, Parkinson's disease, Gaucher disease, schizophrenia and Usher syndrome. Aberrations in chromosome 1 are found in a variety of cancers,



including head and neck cancer, malignant melanoma and multiple myeloma.

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

