



NK-2R Monoclonal Antibody

Catalog No	YP-mAb-12770
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TACR2
Immunogen	The antiserum was produced against synthesized peptide derived from human NK2R. AA range:301-350
Specificity	NK-2R Monoclonal Antibody detects endogenous levels of NK-2R protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Neurokinin A receptor;NK 2 receptor;NK 2R;NK-2 receptor;NK2R;NK-2R;TACR2; Substance-K receptor;NKNAR;TAC2R;SKR;Tachykinin receptor 2
Observed Band	48kD
Calculated Molecular Weight	48kD
Cell Pathway	Cell membrane; Multi-pass membrane protein.
Tissue Specificity	Bone marrow,Lung,Stomach,Trachea,
Function	function:This is a receptor for the tachykinin neuropeptide substance K (neurokinin A). It is associated with G proteins that activate a phosphatidylinositol-calcium second messenger system.,miscellaneous:The rank order of affinity of this receptor to tachykinins is: substance K > neuromedin-K > substance P.,similarity:Belongs to the G-protein coupled receptor 1 family.,
Background	This gene belongs to a family of genes that function as receptors for tachykinins. Receptor affinities are specified by variations in the 5'-end of the sequence. The receptors belonging to this family are characterized by interactions with G proteins and 7 hydrophobic transmembrane regions. This gene encodes the receptor for the tachykinin neuropeptide substance K, also



referred to as neurokinin A. [provided by RefSeq, Jul 2008],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

