



Cleaved-ITI-H1 (D672) Monoclonal Antibody

Catalog No	YP-mAb-15815
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ITIH1
Immunogen	The antiserum was produced against synthesized peptide derived from human ITIH1. AA range:623-672
Specificity	Cleaved-ITI-H1 (D672) Monoclonal Antibody detects endogenous levels of fragment of activated ITI-H1 protein resulting from cleavage adjacent to D672.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Inter-alpha-inhibitor heavy chain 1;Inter-alpha-trypsin inhibitor complex component III;Inter-alpha-trypsin inhibitor heavy chain H1;Serum-derived hyaluronan-associated protein;ITI-H1;ITIH1;IGHEP1;ITI heavy chain H1;ITI-HC1; SHAP
Observed Band	71kD
Calculated Molecular Weight	71kD
Cell Pathway	Secreted.
Tissue Specificity	Blood,Epithelium,Liver,Plasma,Serum,Thymus,
Function	function:Contains a potential peptide which could stimulate a broad spectrum of phagocytotic cells.,function:May act as a carrier of hyaluronan in serum or as a binding protein between hyaluronan and other matrix protein, including those on cell surfaces in tissues to regulate the localization, synthesis and degradation of hyaluronan which are essential to cells undergoing biological processes.,mass spectrometry: PubMed:9677337,polymorphism:There are 3 common alleles; ITIH1*1 with Glu-585/Gln-595, ITIH1*2 with Val-585/Arg-595 and ITIH1*3 with Glu-585/Arg-595.,PTM:Heavy chains are linked to bikunin via chondroitin 4-



sulfate esterified to the alpha-carboxyl of the C-terminal aspartate after propeptide cleavage.,PTM:The S-linked glycan is composed of two units of either Glc or Gal.,similarity:Belongs to the ITIH family.,similarity:Contains 1 VWFA domain.,subunit:I-alpha-I plasma protease in

Background

This gene encodes a member of the inter-alpha-trypsin inhibitor family of proteins. Alternative splicing results in multiple transcript variants, at least one of which encodes a preproprotein that is proteolytically processed to generate the heavy chain of the inter-alpha-trypsin inhibitor complex, which is secreted by hepatocytes into the blood. The heavy chain also interacts with hyaluronan, and this interaction may play a role in ovulation and fertilization, and has been implicated in multiple inflammatory diseases. This gene is present in a gene cluster on chromosome 3. [provided by RefSeq, Nov 2015],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

