



GluR- δ 1 Monoclonal Antibody

Catalog No	YP-mAb-16426
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GRID1
Immunogen	The antiserum was produced against synthesized peptide derived from human GRID1. AA range:831-880
Specificity	GluR- δ 1 Monoclonal Antibody detects endogenous levels of GluR- δ 1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	$\geq 90\%$
Storage Stability	-20°C/1 year
Synonyms	delta-2;GluD2;GLUR D2;GluR delta 2;GluR delta-2 subunit;GLURD 2;GLURD2; Glutamate receptor delta 2 subunit;Glutamate receptor ionotropic;Glutamate receptor ionotropic delta 2;GRID 2;Grid2;GRID2_HUMAN;MGC11722; MGC11723;MGC11724;GluR- δ 1;GRID1;Glutamate receptor delta-1 subunit; KIAA1220;GluR delta-1 subunit
Observed Band	95kD
Calculated Molecular Weight	95kD
Cell Pathway	Cell membrane ; Multi-pass membrane protein . Cell junction, synapse, postsynaptic cell membrane ; Multi-pass membrane protein .
Tissue Specificity	Brain,Uterus,
Function	function:Receptor for glutamate. L-glutamate acts as an excitatory neurotransmitter at many synapses in the central nervous system. The postsynaptic actions of Glu are mediated by a variety of receptors that are named according to their selective agonists.,similarity:Belongs to the glutamate-gated ion channel (TC 1.A.10) family.,



Background

This gene encodes a subunit of glutamate receptor channels. These channels mediate most of the fast excitatory synaptic transmission in the central nervous system and play key roles in synaptic plasticity.[provided by RefSeq, Jan 2009],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using GluR- δ 1 Monoclonal Antibody