



Caspase-10 B/C Monoclonal Antibody

Catalog No	YP-mAb-00342
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CASP10
Immunogen	The antiserum was produced against synthesized peptide derived from human Caspase 10. AA range:430-479
Specificity	Caspase-10 B/C Monoclonal Antibody detects endogenous levels of Caspase-10 B/C protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CASP10;Caspase10;ALPS2;apoptotic protease Mch-4;CASP-10;CASPA; caspase-10;FADD-like ICE2;FAS-associated death domain protein interleukin-1B-convert
Observed Band	60
Calculated Molecular Weight	58994
Cell Pathway	cytosol,CD95 death-inducing signaling complex,ripiptosome,
Tissue Specificity	Detectable in most tissues. Lowest expression is seen in brain, kidney, prostate, testis and colon.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images