



POLR1E Monoclonal Antibody

Catalog No	YP-mAb-01946
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	POLR1E
Immunogen	Synthesized peptide derived from the Internal region of human POLR1E.
Specificity	POLR1E Monoclonal Antibody detects endogenous levels of POLR1E protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	A49;DNA-directed RNA polymerase I subunit E;DNA-directed RNA polymerase I subunit RPA49;PAF53;POLR;POLR1E;PRAF1;RNA polymerase I subunit A49; RNA polymerase I-associated factor 1;RNA polymerase I-associated factor 53
Observed Band	49kD
Calculated Molecular Weight	49kD
Cell Pathway	Nucleus, nucleolus .
Tissue Specificity	Cervix,Lung,Placenta,
Function	function:DNA-dependent RNA polymerase catalyzes the transcription of DNA into RNA using the four ribonucleoside triphosphates as substrates. Component of RNA polymerase I which synthesizes ribosomal RNA precursors. Appears to be involved in the formation of the initiation complex at the promoter by mediating the interaction between Pol I and UBTF/UBF.,similarity:Belongs to the eukaryotic RPA49/POLR1E RNA polymerase subunit family., subunit:Component of the RNA polymerase I (Pol I) complex consisting of at least 13 subunits (By similarity). Interacts with PAF49/CD3EAP. Also binds UBTF/UBF.,

**Background**

function:DNA-dependent RNA polymerase catalyzes the transcription of DNA into RNA using the four ribonucleoside triphosphates as substrates. Component of RNA polymerase I which synthesizes ribosomal RNA precursors. Appears to be involved in the formation of the initiation complex at the promoter by mediating the interaction between Pol I and UBTF/UBF.,similarity:Belongs to the eukaryotic RPA49/POLR1E RNA polymerase subunit family., subunit:Component of the RNA polymerase I (Pol I) complex consisting of at least 13 subunits (By similarity). Interacts with PAF49/CD3EAP. Also binds UBTF/UBF.,

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images