



ZNF446 Monoclonal Antibody

Catalog No	YP-mAb-02189
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZNF446
Immunogen	The antiserum was produced against synthesized peptide derived from human ZNF446. AA range:141-190
Specificity	ZNF446 Monoclonal Antibody detects endogenous levels of ZNF446 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZSCAN30;ZKSCAN20;ZNF446;Zinc finger protein 446;Zinc finger protein with KRAB and SCAN domains 20
Observed Band	49-55 kDa
Calculated Molecular Weight	422 aa, 46 kDa
Cell Pathway	Nucleus .
Tissue Specificity	Placenta,Testis,
Function	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 1 SCAN box domain.,similarity:Contains 3 C2H2-type zinc fingers.,
Background	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 1 SCAN box domain.,similarity:Contains 3 C2H2-type zinc fingers.,

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images