



DHRS11 Monoclonal Antibody

Catalog No	YP-mAb-02618
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DHRS11
Immunogen	The antiserum was produced against synthesized peptide derived from human DHRS11. AA range:111-160
Specificity	DHRS11 Monoclonal Antibody detects endogenous levels of DHRS11 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ARPG836;DHRS11;SDR24C1;Dehydrogenase/reductase SDR family member 11
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	Isoform 1: Ubiquitously expressed, with highest levels in testis, small intestine, colon, kidney, brain and heart. Isoform 3: Expressed in brain, heart and skeletal muscle.
Function	similarity:Belongs to the short-chain dehydrogenases/reductases (SDR) family., subunit:Homotetramer.,
Background	similarity:Belongs to the short-chain dehydrogenases/reductases (SDR) family., subunit:Homotetramer.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images