



PDHK1 Monoclonal Antibody

Catalog No	YP-mAb-02742
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PDPK1
Immunogen	The antiserum was produced against synthesized peptide derived from human PDK1. AA range:1-50
Specificity	PDK1 Monoclonal Antibody detects endogenous levels of PDK1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PDPK1;3-phosphoinositide-dependent protein kinase 1;PDK1;hPDK1
Observed Band	60
Calculated Molecular Weight	63152
Cell Pathway	Cytoplasm. Nucleus. Cell membrane; Peripheral membrane protein. Cell junction, focal adhesion. Tyrosine phosphorylation seems to occur only at the cell membrane. Translocates to the cell membrane following insulin stimulation by a mechanism that involves binding to GRB14 and INSR. SRC and HSP90 promote its localization to the cell membrane. Its nuclear localization is dependent on its association with PTPN6 and its phosphorylation at Ser-396. Restricted to the nucleus in neuronal cells while in non-neuronal cells it is found in the cytoplasm. The Ser-241 phosphorylated form is distributed along the perinuclear region in neuronal cells while in non-neuronal cells it is found in both the nucleus and the cytoplasm. IGF1 transiently increases phosphorylation at Ser-241 of neuronal PDPK1, resulting in its translocation to other cellular compartments. The tyrosine-phosphorylated form colocalizes with PTK2B in focal adhesions after angiotensin II stimulation.



Tissue Specificity	Appears to be expressed ubiquitously. The Tyr-9 phosphorylated form is markedly increased in diseased tissue compared with normal tissue from lung, liver, colon and breast.
---------------------------	--

Function	
-----------------	--

Background	
-------------------	--

matters needing attention	Avoid repeated freezing and thawing!
----------------------------------	--------------------------------------

Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.
--------------------------	---

Products Images