



USP32 Monoclonal Antibody

Catalog No	YP-mAb-02833
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	USP32
Immunogen	The antiserum was produced against synthesized peptide derived from human USP32. AA range:271-320
Specificity	USP32 Monoclonal Antibody detects endogenous levels of USP32 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	USP10;Ubiquitin-specific-processing protease 32;Ubiquitin carboxyl-terminal hydrolase 32;Renal carcinoma antigen NY-REN-60;NY REN 60;USP32;Deubiquitinating enzyme 32;Ubiquitin thioesterase 32
Observed Band	180kD
Calculated Molecular Weight	180kD
Cell Pathway	Membrane ; Lipid-anchor . Golgi apparatus .
Tissue Specificity	Brain,Epithelium,Peripheral Nervous System,Skin,Testis,Uterus,
Function	catalytic activity:Ubiquitin C-terminal thioester + H(2)O = ubiquitin + a thiol., similarity:Belongs to the peptidase C19 family.,similarity:Contains 1 DUSP domain.,similarity:Contains 3 EF-hand domains.,
Background	catalytic activity:Ubiquitin C-terminal thioester + H(2)O = ubiquitin + a thiol., similarity:Belongs to the peptidase C19 family.,similarity:Contains 1 DUSP domain.,similarity:Contains 3 EF-hand domains.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images