



CLIP-170 Monoclonal Antibody

Catalog No	YP-mAb-03102
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CLIP1
Immunogen	The antiserum was produced against synthesized peptide derived from human CLIP1. AA range:1291-1340
Specificity	CLIP-170 Monoclonal Antibody detects endogenous levels of CLIP-170 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CAP-Gly domain-containing linker protein 1;CLIP-170;CYLN1;Cytoplasmic linker protein 170 alpha-2;Reed-Sternberg intermediate filament-associated protein; CLIP1;RSN;Cytoplasmic linker protein 1;Restin
Observed Band	161kD
Calculated Molecular Weight	161kD
Cell Pathway	Cytoplasm . Cytoplasm, cytoskeleton . Cytoplasmic vesicle membrane ; Peripheral membrane protein; Cytoplasmic side. Cell projection, ruffle . Localizes to microtubule plus ends (PubMed:21646404, PubMed:17889670). Localizes preferentially to the ends of tyrosinated microtubules (By similarity). Accumulates in plasma membrane regions with ruffling and protrusions. Associates with the membranes of intermediate macropinocytic vesicles (PubMed:12433698). .
Tissue Specificity	Detected in dendritic cells (at protein level). Highly expressed in the Reed-Sternberg cells of Hodgkin disease.
Function	function:Seems to be a intermediate filament associated protein that links endocytic vesicles to microtubules.,similarity:Contains 2 CAP-Gly domains., subcellular location:Associated with the cytoskeleton.,tissue specificity:Highly expressed in the Reed-Sternberg cells of Hodgkin's disease.,

**Background**

The protein encoded by this gene links endocytic vesicles to microtubules. This gene is highly expressed in Reed-Sternberg cells of Hodgkin disease. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2011],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images