



KIF13B Monoclonal Antibody

Catalog No	YP-mAb-03143
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KIF13B
Immunogen	The antiserum was produced against synthesized peptide derived from human KIF13B. AA range:111-160
Specificity	KIF13B Monoclonal Antibody detects endogenous levels of KIF13B protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GAKIN;KIAA0639;kinesin family member 13B;KIF13B;Kinesin-like protein KIF13B;Kinesin-like protein GAKIN
Observed Band	202kD
Calculated Molecular Weight	202kD
Cell Pathway	Cytoplasm, cytoskeleton . Cell projection, axon . accumulates at the distal part of the microtubules in the tips of axons, but not of dendrites.
Tissue Specificity	Ubiquitous.
Function	function:May be involved in reorganization of the cortical cytoskeleton. May be functionally important for the intracellular trafficking of MAGUKs and associated protein complexes.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the kinesin-like protein family.,similarity:Contains 1 CAP-Gly domain.,similarity:Contains 1 FHA domain.,similarity:Contains 1 kinesin-motor domain.,subunit:Binds to DLG1 and DLG4.,tissue specificity:Ubiquitous.,

**Background**

function: May be involved in reorganization of the cortical cytoskeleton. May be functionally important for the intracellular trafficking of MAGUKs and associated protein complexes., PTM: Phosphorylated upon DNA damage, probably by ATM or ATR., similarity: Belongs to the kinesin-like protein family., similarity: Contains 1 CAP-Gly domain., similarity: Contains 1 FHA domain., similarity: Contains 1 kinesin-motor domain., subunit: Binds to DLG1 and DLG4., tissue specificity: Ubiquitous.,

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images