



Myosin VA Monoclonal Antibody

Catalog No	YP-mAb-03166
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MYO5A
Immunogen	The antiserum was produced against synthesized peptide derived from human MYO5A. AA range:1784-1833
Specificity	Myosin VA Monoclonal Antibody detects endogenous levels of Myosin VA protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GS1;MYH12;MYO5A;Myosin 12;Myosin heavy chain 12;Myosin Va;Myoxin;MYR12;Unconventional myosin-Va;Dilute myosin heavy chain;non-muscle;Myosin-12
Observed Band	220kD
Calculated Molecular Weight	220kD
Cell Pathway	ruffle,photoreceptor outer segment,cytoplasm,lysosome,early endosome,late endosome,peroxisome,endoplasmic reticulum,Golgi apparatus,cytosol,intermediate filament,actin filament,membrane,myosin complex,gr
Tissue Specificity	Detected in melanocytes.
Function	disease:Defects in MYO5A are a cause of Elejalde syndrome [MIM:256710]; also known as neuroectodermal melanolysosomal disease. Elejalde syndrome is an autosomal recessive condition characterized by skin hypopigmentation, the presence of large clumps of pigment in hair shafts, silvery-gray hair, accumulation of melanosomes in melanocytes and primary neurological abnormalities. Elejalde syndrome may be the same entity as Griscelli syndrome type I.,disease:Defects in MYO5A are a cause of Griscelli syndrome type-1 (GS1) [MIM:214450]; also known as Griscelli syndrome with primary neurologic



impairment. Griscelli syndrome is a rare autosomal recessive disorder that results in pigmentary dilution of the skin and hair, the presence of large clumps of pigment in hair shafts, silvery-gray hair and accumulation of melanosomes in melanocytes. GS1 patients show developmental delay, hypotonia and ment

Background

This gene is one of three myosin V heavy-chain genes, belonging to the myosin gene superfamily. Myosin V is a class of actin-based motor proteins involved in cytoplasmic vesicle transport and anchorage, spindle-pole alignment and mRNA translocation. The protein encoded by this gene is abundant in melanocytes and nerve cells. Mutations in this gene cause Griscelli syndrome type-1 (GS1), Griscelli syndrome type-3 (GS3) and neuroectodermal melanolyosomal disease, or Elejalde disease. Multiple alternatively spliced transcript variants encoding different isoforms have been reported, but the full-length nature of some variants has not been determined. [provided by RefSeq, Dec 2008],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images