



Myosin XVI Monoclonal Antibody

Catalog No	YP-mAb-03167
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MYO16
Immunogen	Synthesized peptide derived from Myosin XVI . at AA range: 1050-1130
Specificity	Myosin XVI Monoclonal Antibody detects endogenous levels of Myosin XVI protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIAA0865;MYO16;Myo16b;myosin XVI;MYR8;Unconventional myosin 16; Unconventional myosin-XVI;NYAP3;Neuronal tyrosine-phosphorylated phosphoinositide-3-kinase adapter 3;Unconventional myosin-16
Observed Band	200kD
Calculated Molecular Weight	200kD
Cell Pathway	Cytoplasm . Found in puncta in soma and processes of astrocytes and dissociated cerebellar cells with the morphology of migrating granule cells. .
Tissue Specificity	Amygdala,Brain,
Function	function:Myosins are actin-based motor molecules with ATPase activity. Unconventional myosins serve in intracellular movements. Their highly divergent tails are presumed to bind to membranous compartments, which would be moved relative to actin filaments. May be involved in targeting of the catalytic subunit of protein phosphatase 1 during brain development.,similarity:Contains 1 IQ domain.,similarity:Contains 2 myosin head-like domains.,similarity:Contains 7 ANK repeats.,subcellular location:Found in puncta in soma and processes of astrocytes and dissociated cerebellar cells with the morphology of migrating granule cells.,subunit:Binds PPP1CA and/or PPP1CC. Binds F-actin in an ATP-



sensitive manner.,

Background

function:Myosins are actin-based motor molecules with ATPase activity. Unconventional myosins serve in intracellular movements. Their highly divergent tails are presumed to bind to membranous compartments, which would be moved relative to actin filaments. May be involved in targeting of the catalytic subunit of protein phosphatase 1 during brain development.,similarity:Contains 1 IQ domain.,similarity:Contains 2 myosin head-like domains.,similarity:Contains 7 ANK repeats.,subcellular location:Found in puncta in soma and processes of astrocytes and dissociated cerebellar cells with the morphology of migrating granule cells.,subunit:Binds PPP1CA and/or PPP1CC. Binds F-actin in an ATP-sensitive manner.,

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images