



Gab 2/4 Monoclonal Antibody

Catalog No	YP-mAb-03885
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GAB2/GAB4
Immunogen	The antiserum was produced against synthesized peptide derived from human GAB2. AA range:21-70
Specificity	Gab 2/4 Monoclonal Antibody detects endogenous levels of Gab 2/4 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Gab 2/4;GAB2/GAB4;GRB2-associated-binding protein 2/4;GAB4;GRB2-associated-binding protein 4;GRB2-associated binder 2-like;GAB2-like;GRB2-associated binder 4;GRB2-associated-binding protein 2-like;Growth factor receptor bound protein 2-associated protein 4;GAB2;KIAA0571;GRB2-associate
Observed Band	62kD
Calculated Molecular Weight	62kD
Cell Pathway	
Tissue Specificity	
Function	similarity:Belongs to the GAB family.,similarity:Contains 1 PH domain.,
Background	similarity:Belongs to the GAB family.,similarity:Contains 1 PH domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images