



LONP2 Monoclonal Antibody

Catalog No	YP-mAb-03946
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LONP2
Immunogen	The antiserum was produced against synthesized peptide derived from human LONP2. AA range:780-829
Specificity	LONP2 Monoclonal Antibody detects endogenous levels of LONP2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EC:3.4.21.53;Lon protease 2;Lon protease homolog 2;peroxisomal;Lon protease-like protein 2;LONP;LONP2;Lon protease homolog 2 peroxisomal; Peroxisomal Lon protease
Observed Band	95kD
Calculated Molecular Weight	95kD
Cell Pathway	Peroxisome matrix .
Tissue Specificity	Widely expressed, with high levels in the liver, kidney and pancreas.
Function	function:Probable serine protease which may be involved in peroxisome biogenesis.,similarity:Belongs to the peptidase S16 family.,
Background	function:Probable serine protease which may be involved in peroxisome biogenesis.,similarity:Belongs to the peptidase S16 family.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images