



SH3RF3 Monoclonal Antibody

Catalog No	YP-mAb-04207
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SH3RF3
Immunogen	Synthesized peptide derived from the C-terminal region of human SH3RF3.
Specificity	SH3RF3 Monoclonal Antibody detects endogenous levels of SH3RF3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	POSH2;SH3MD4;SH3RF3;SH3 domain-containing RING finger protein 3;Plenty of SH3s 2;SH3 multiple domains protein 4
Observed Band	93 kDa
Calculated Molecular Weight	92kD
Cell Pathway	
Tissue Specificity	Brain,Spleen,
Function	domain:The third SH3 domain mediates interaction with PAK2., similarity:Belongs to the SH3RF family.,similarity:Contains 1 RING-type zinc finger.,similarity:Contains 4 SH3 domains.,subunit:Interacts with PAK2.,
Background	domain:The third SH3 domain mediates interaction with PAK2., similarity:Belongs to the SH3RF family.,similarity:Contains 1 RING-type zinc finger.,similarity:Contains 4 SH3 domains.,subunit:Interacts with PAK2.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images