



BCAR1 mouse mAb

Catalog No	YP-mAb-04363
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	BCAR1 CAS CASS1 CRKAS
Immunogen	Synthesized peptide derived from human BCAR1
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat BCAR1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BCAR1 CAS CASS1 CRKAS;BCAR1;Breast cancer anti-estrogen resistance protein 1 (CRK-associated substrate;Cas scaffolding protein family member 1; p130cas)
Observed Band	91868DA
Calculated Molecular Weight	95700
Cell Pathway	Cell junction, focal adhesion . Cytoplasm . Cell projection, axon . Unphosphorylated form localizes in the cytoplasm (By similarity). Localizes to focal adhesion sites following integrin engagement (By similarity). .
Tissue Specificity	Widely expressed with an abundant expression in the testis. Low level of expression seen in the liver, thymus, and peripheral blood leukocytes. The protein has been detected in a B-cell line.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images