



CATF (Cleaved-Ala271) mouse mAb

Catalog No	YP-mAb-04366
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CTSF
Immunogen	Synthesized peptide derived from human CATF (Cleaved-Ala271)
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat CATF (Cleaved-Ala271, protein was cleaved amino acid sequence between 270-271)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CATF;CTSF;CATF (Cleaved-Ala271);Cathepsin F (CATSF;EC 3.4.22.41)
Observed Band	23 55kD
Calculated Molecular Weight	23 55kD
Cell Pathway	Lysosome.
Tissue Specificity	High expression levels in heart, skeletal muscle, brain, testis and ovary; moderate levels in prostate, placenta, liver and colon; and no detectable expression in peripheral leukocytes and thymus.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images