



CATW (Cleaved-Val128) mouse mAb

Catalog No	YP-mAb-04369
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CTSW
Immunogen	Synthesized peptide derived from human CATW (Cleaved-Val128)
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat CATW (Cleaved-Val128, protein was cleaved amino acid sequence between 127-128)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CATW;CTSW;CATW (Cleaved-Val128);Cathepsin W (EC 3.4.22.-;Lymphopain)
Observed Band	28 42kD
Calculated Molecular Weight	28 42kD
Cell Pathway	Endoplasmic reticulum .
Tissue Specificity	Expressed predominantly in natural killer cells, and in cytotoxic T cells.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images