



Collagen XV α 1 (Cleaved-Asn1212) mouse mAb

Catalog No	YP-mAb-04380
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	COL15A1
Immunogen	Synthesized peptide derived from human Collagen XV α 1 (Cleaved-Asn1212)
Specificity	This antibody detects endogenous levels of Human,Mouse Collagen XV α 1 (Cleaved-Asn1212, protein was cleaved amino acid sequence between 1211-1212)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	$\geq 90\%$
Storage Stability	-20°C/1 year
Synonyms	Collagen XV α 1;COL15A1;Collagen XV α 1 (Cleaved-Asn1212);Collagen alpha-1(XV) chain [Cleaved into: Endostatin (Endostatin-XV;Related to endostatin; Restin)]
Observed Band	20 150kD
Calculated Molecular Weight	20 150kD
Cell Pathway	Secreted, extracellular space, extracellular matrix .
Tissue Specificity	Expressed predominantly in internal organs such as adrenal gland, pancreas and kidney.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images